

다양한 형광 단백질 검출을 위한

Living Colors® Antibodies

AcGFP1, EGFP, GFP, ZsGreen1



형광단백질	검출가능한 항체 제품명	Code	용도*
AcGFP1	A.v. Monoclonal Antibody (JL-8)	632381	WB, IP, IC
	Full Length GFP Polyclonal Antibody	632592	WB, IP, IC
EGFP**	EGFP Monoclonal Antibody	632569	WB, IP, IC
	Full-Length GFP Polyclonal Antibody	632592	WB, IP, IC
GFP & GFPuv	Full Length GFP Polyclonal Antibody	632592	WB, IL
ZsGreen1	Full-Length ZsGreen Monoclonal Antibody	632598	WB
	Full-Length ZsGreen1 Polyclonal Antibody	632474	WB, IP, IL

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry, IL = immunolabeling

** Not sold by Takara Bio

AmCyan1, ECFP



형광단백질	검출가능한 항체 제품명	Code	용도*
AmCyan1	Anti-RCFP Polyclonal Pan Antibody	632475	WB
ECFP**	A.v. Monoclonal Antibody (JL-8)	632381	WB, IP, IC
	Full-Length GFP Polyclonal Antibody***	632592	WB, IP, IC

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry

** Not sold by Takara Bio

*** ECFP is highly related to EGFP and is recognized by the same antibodies

mBanana, EYFP, ZsYellow



형광단백질	검출가능한 항체 제품명	Code	용도*
mBanana	DsRed Polyclonal Antibody	632496	WB
EYFP**	A.v. Monoclonal Antibody (JL-8)	632381	WB, IP, IC
	Full-Length GFP Polyclonal Antibody***	632592	WB, IP, IC
ZsYellow1	Anti-RCFP Polyclonal Pan Antibody	632475	WB, IL

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry, IL = immunolabeling

** Not sold by Takara Bio

*** EYFP is highly related to EGFP so is detected by the same antibodies

mOrange2



형광단백질	검출가능한 항체 제품명	Code	용도*
mOrange2	mCherry Monoclonal Antibody	632543	WB, IC

* WB = Western Blot, IC = immunocytochemistry


AsRed, mCherry, DsRed, tdTomato



형광단백질	검출가능한 항체 제품명	Code	용도*
AsRed2	Anti-RCFP Polyclonal Pan Antibody	632475	WB, IL
mCherry	mCherry Monoclonal Antibody	632543	WB, IP, IL
	DsRed Polyclonal Antibody	632496	WB, IP, IL
Photoactivatable mCherry	mCherry Monoclonal Antibody	632543	WB, IP, IL
DsRed2	DsRed Monoclonal Antibody	632392	WB
	DsRed Polyclonal Antibody	632496	WB, IP, IL
DsRed-Express	DsRed Monoclonal Antibody	632392	WB
	DsRed Polyclonal Antibody	632496	WB, IP, IL
DsRed-Express2	DsRed Monoclonal Antibody	632392	WB
	DsRed Polyclonal Antibody	632496	WB, IP, IL
DsRed-Monomer	Anti-RCFP Polyclonal Pan Antibody	632475	WB, IL
mStrawberry	mCherry Monoclonal Antibody	632543	WB, IP, IL
	DsRed Polyclonal Antibody	632496	WB, IP, IL
tdTomato	DsRed Polyclonal Antibody	632496	WB, IP, IL
Timer (DsRed1-E5)	DsRed Polyclonal Antibody	632496	WB, IP, IL

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry, IL = immunolabeling

HcRed1, E2-Crimson, mPlum, mRaspberry



형광단백질	검출가능한 항체 제품명	Code	용도*
HcRed1	Anti-RCFP Polyclonal Pan Antibody	632475	WB
E2-Crimson	DsRed Monoclonal Antibody	632392	WB
	DsRed Polyclonal Antibody	632496	WB
mPlum	mCherry Monoclonal Antibody	632543	WB, IC
	DsRed Polyclonal Antibody	632496	WB
mRaspberry	DsRed Monoclonal Antibody	632392	WB
	DsRed Polyclonal Antibody	632496	WB

* WB = Western Blot, IC = immunocytochemistry

같이 쓰면 좋은 제품

형광벡터 클로닝 효율 100%를 위한,

The Best Cloning Kit

In-Fusion® Snap Assembly Master Mix (Code 638947 외)

- The first product of seamless, directional cloning

