## **CHOgro® Complex Formation Solution (100 ml)**

## **Product Data Sheet**

Product information for MIR 6210
SDS and Certificate of Analysis available at mirusbio.com/literature



## **SPECIFICATIONS**

Storage	Store CHOgro® Complex Formation Solution at 4°C, protected from light.
Product Guarantee	As indicated on the product label, when properly stored and handled.

Product Description: CHOgro® Complex Formation Solution is a chemically defined, hydrolysate-

free and animal-origin-free medium. It is designed for use with *Trans*IT-PRO® Transfection Reagent within the CHOgro® Expression System and CHOgro® High Yield Expression System for high and reproducible protein yields in suspension CHO cells. CHOgro® Complex Formation Solution is supplied in a 100 ml bottle format and is also a component of the complete CHOgro® Expression System (MIR 6260) and CHOgro® High Yield Expression System (MIR 6270).

System (min sess) and enoges mgn near expression system (min sers).

Use the CHOgro® Complex Formation Solution to facilitate complex formation between DNA and transfection reagent prior to transfection. See the CHOgro®

High Yield Expression System Full Protocol for a full protocol.

Product Size: 100 ml

**Product Usage:** 

Form/Appearance: Clear, yellow solution.

Precautions: Review the product-specific Safety Data Sheet (SDS) for additional information

before handling this product.

For Research Use Only



SDS and Certificate of Analysis available at <u>mirusbio.com/literature</u>

©1996-2024 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners. For terms and conditions, visit www.mirusbio.com

Rev1 040CT2024